

The neria™ infusion set line

Product Specifications	neria™	neria™ detach	neria™ multi	neria™ soft	neria™ soft 90
Cannula	Stainless Steel	Stainless Steel	Stainless Steel	Soft	Soft
Insertion angle	90 degrees	90 degrees	90 degrees	20-45 degrees	90 degrees
Disconnection feature	No	Yes	No	Yes at site	Yes at site
Needle gauge	G27 or G29	G27	G27	G27	G27
Needle lengths	6, 8 or 10 mm	6, 8 or 10 mm	Bi-furcated: 8, 10 mm Tri-furcated: 8, 10, 12 mm Quad-furcated: 8, 10, 12 mm	13 or 17 mm	6 or 9 mm
Tubing lengths	60 cm, 80 cm, 110 cm	60 cm, 80 cm, 110 cm	Total 90 cm	30, 60, 80, 110 cm	30, 60, 110 cm
Adhesive	Built-in	Built-in	Built-in	Built-in	Built-in
Priming volumes	60 cm - 0.10 ml 80 cm - 0.13 ml 110 cm - 0.18 ml	60 cm - 0.10 ml 80 cm - 0.13 ml 110 cm - 0.18 ml	Bi-furcated: 0.36 ml Tri-furcated: 0.43 ml Quad-furcated: 0.51 ml	30 cm - 0.05 ml 60 cm - 0.10 ml 80 cm - 0.13 ml 110 cm - 0.18 ml	30 cm - 0.05 ml 60 cm - 0.10 ml 110 cm - 0.18 ml