

Product Specifications	neria™	neria™ detach	neria™ multi	neria™ soft	neria™ soft 90	neria™ guard
Cannula	Stainless Steel	Stainless Steel	Stainless Steel	Soft	Soft	Soft
Insertion angle	90 degrees	90 degrees	90 degrees	20-45 degrees	90 degrees	90 degrees
Disconnection feature	No	Yes	No	Yes at site	Yes at site	Yes at site
Needle/introducer needle gauge	G27 or G29	G27	G27	G27	G27	G27
Needle lengths	6, 8 or 10 mm	6, 8 or 10 mm	Bi-furcated: 8, 10 mm Tri-furcated: 8, 10, 12 mm Quad-furcated: 8, 10, 12 mm	13 or 17 mm	6 or 9 mm	6 or 9 mm
Tubing lengths	60, 80, 110 cm	60, 80, 110 cm	Total 90 cm	30, 60, 80, 110 cm	30, 60, 110 cm	12, 60, 80, 110 cm
Adhesive	Built-in	Built-in	Built-in	Built-in	Built-in	Built-in
Needle safety	-	-	-	Yes	Yes	Yes
Priming volumes	60 cm ~ 0.10 ml 80 cm ~ 0.12 ml 110 cm ~ 0.15 ml	60 cm ~ 0.10 ml 80 cm ~ 0.12 ml 110 cm ~ 0.15 ml	Bi-furcated ~ 0.36 ml Tri-furcated ~ 0.43 ml Quad-furcated ~ 0.51 ml	30 cm ~ 0.06 ml 60 cm ~ 0.10 ml 80 cm ~ 0.12 ml 110 cm ~ 0.15 ml	30 cm ~ 0.06 ml 60 cm ~ 0.10 ml 110 cm ~ 0.15 ml	12 cm ~ 0.04 ml 60 cm ~ 0.10 ml 80 cm ~ 0.12 ml 110 cm ~ 0.15 ml